



KLINIPATH Microscope Slides

Description:	SLIDES, Ground Edges ,Clipped Corner Slide； Double frosted, etched both sides one end； Frosted Area: 19.5 x 25 (±0.5)mm;
Purpose:	Used for hematology ,microbiology,histology pathology
Materials:	Standard Glass; thinckness 1.0 (-0.05,+0.1)mm
Appearance:	Pre-cleaned,clear;no marks,prints or contamination,etc
Features:	The smooth edges provide the safety of the operation
	For marking with pencil ,mark pens and labeling system
Temperature	Normal laboratory temperature and humidity
Uses:	For human body blood, body fluids test, and paraffin section
Shelf life：	60 Months

KLINIPATH Microscope Slides

EU Part No.	Dimension				Frosted or Painted	Paint Color	Corner Angle	Edges	Packaging	QTY/CASE
	A (mm)	B (mm)	C (mm)	D (mm)						
			-0.05/+0.1	(+/-0.5)						
KP-3004	75.7~76.0	25.0~25.3	1	N/A	Double frosted ends	N/A	45° (small corner)	Smooth	50pcs / Box, 50 Boxes / Case	2500