n-Heptane 99%
ChromAR HPLC

Material No: 5410-25
Effective from: 2014/03/13

Product Specification

Test	Specification
Assay (GC)	>= 99 %
Water (KF)	<= 0.01 %
Residue after Evaporation	<= 0.0005 %
Titrable Acid (meq/g)	<= 0.0005
Ultraviolet Absorbance (1.00-cm cell vs water)	
at 210 nm	<= 0.6
at 220 nm	<= 0.2
at 230 nm	<= 0.05
at 240 nm	<= 0.02
at 250 nm	<= 0.006

For Laboratory, Research or Manufacturing use.
Suitable for Liquid Chromatography and UV-Spectrophotometry

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000