Methanol
HPLC
Gradient Grade: For use in Liquid Chromatography
(HPLC & UHPLC) & Spectrophotometry
(methyl alcohol)

Material No.: 9093-02
:
Revision No.: 0

Product Specification

Test	Specification
Assay (CH_3OH) (by GC, corrected for wate	≥ 99.9 %
Ultraviolet Absorbance (1.00-cm cell v	≤ 0.01
Ultraviolet Absorbance (1.00-cm cell v	≤ 0.15
Ultraviolet Absorbance (1.00-cm cell v	≤ 205
Gradient Elution Test (a.u.) - 254 nm	≤ 0.002
Fluorescence Trace Impurities, measured	≤ 0.3 ppb
Fluorescence Trace Impurities, measured	≤ 1.0 ppb
Acetone	≤ 0.001 %
Residue after Evaporation	≤ 1.0 ppm
Titrable Acid (µeq/g)	≤ 0.30
Titrable Base (µeq/g)	≤ 0.10
Water (by KF, coulometric)	≤ 0.05 %

For Laboratory,Research,or Manufacturing Use
Filtered through a 0.2 micron filter.