# Laboratory bowls, flat form



Laboratory bowl, flat form

| Material | C (Height) | D/Ø | l (Liter) | Article no. |
|---|---|---|---|---|
| 18/10 Stainless | 75 | 165 | 1 | 8590 |
| 18/10 Stainless | 75 | 220 | 2 | 8591 |
| 18/10 Stainless | 95 | 240 | 3 | 8592 |
| 18/10 Stainless | 120 | 260 | 4,5 | 8593 |
| 18/10 Stainless | 120 | 305 | 6 | 8594 |
| 18/10 Stainless | 145 | 320 | 8 | 8595 |
| 18/10 Stainless | 135 | 380 | 11 | 8596 |
| 18/10 Stainless | 160 | 400 | 14 | 8597 |