Spécifications du produit

Produit	94006.1016
Description du produit	Ipasept
Grade / qualité	GAMMA-IRRADIATED 25 kGy 0,2 µm filtrated
Numéro CAS	67-63-0
Formule moléculaire	$CH_3CHOHCH_3$

Analyses	Spécifications
Titre (Propanol-2) (V/V)	68.0 - 72.0 %
Titre (Propanol-2) (V/V) [USP]	Conforme USP
Aspect	Liquide limpide incolore
Acidité	Conforme USP
Matières premières	Conforme Ph.Eur.
Densité (20/4)	0.872 - 0.883
Matières non volatiles	≤ 100 mg/l

Signature

Nous certifions que ce lot est conforme aux spécifications listées ci-dessus.

Ce document a été produit électroniquement et est validé sans signature.

Jessie Pijck, Head of Laboratory - Haasrode
VWR International bv; Geldenaaksebaan 464; BE-3001 Leuven; Belgium

The irradiation process has been validated according to ISO 11137-2:2013, VDMAX 25 method.
Starting form batch 14F a sterility assurance level of 10-6 can be guaranteed on the entire content of the inner bag.

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.