

ThermoFisher SCIENTIFIC

Product Specifications

Product **153230250**
Giemsa Stain,pure,high purity biological stai n

General Product Data

Version	2
CAS No	51811-82-6
Molecular weight	0.00
Molecular formula	
Linear formula	

Product Specifications

Appearance (Color)	Dark green to black
Appearance (Form)	Crystalline powder or crystals
Infrared spectrum	Conforms
Loss on drying	=<10 % (110°C)
UV Lambda max	640 to 650 nm (in MeOH)
Specific absorbance A (1%/1cm)	>=950 (in MeOH) On dried substance
UV Lambda max	520 to 525 nm (in MeOH)
Specific absorbance A (1%/1cm)	>=600 (in MeOH) On dried substance

Product Specifications are subject to amendment and may change over time.

Acros Organics BV

ENA23, zone 1, nr 1350, Janssen Pharmaceuticalaan 3a, B-2440 Geel, Belgium
Tel +32 14/57.52.11 - Fax +32 14/59.34.34 Internet: http://www.thermofisher.com
1 Reagent Lane, Fair Lawn, NJ 07410,USA
Fax 201-796-1329