BAKERBOND spe, Aromatic Sulfonic Acid

BAKERBOND spe

Disposable Extraction Columns

Material No: 7090-01
Effective from: 2016/03/23

# Product Specification

| Test | Specification |
|---|---|
| The following results are based on the bulk bonded phase: | |
| Carbon (C)(%) | For Information Only |
| Sulfur (S)(ug/ml) | For Information Only |
| Surface Coverage (ug/m2), Based on %C | 180 - 245 |
| Surface Coverage (ug/m2), Based on %S | 50 - 80 |
| The following results are based on the silica gel raw material: | |
| Average Particle Diameter, um (APD) | 47 - 60 |
| The following results are based on the solid phase extraction columns: | |
| Extraction Residue | <= 0.2 % |

For Laboratory, Research or Manufacturing use.
1mL Solid Phase Extraction Columns
100 mg per Column

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000