Hydrochloric Acid, 5N Volumetric Solution
BAKER ANALYZED® Reagent

Material No.: 5618-07
:
Revision No.: 0

Product Specification

Test	Specification
Normality	4.975 - 5.025
Sulfate (SO_4)	≤ 1 ppm
Ammonium (NH_4)	≤ 3 ppm
Trace Impurities - Heavy Metals (as Pb)	≤ 0.5 ppm
Trace Impurities - Iron (Fe)	≤ 0.1 ppm
Appearance	Passes Test

For Laboratory,Research,or Manufacturing Use
Standardization at 25°C traceable to NIST Standard Reference Material.

Packaging Site: Paris Mfg Ctr & DC