Technical Data Sheet

Edition Date: 2016-10-04
Revision Date : 2018-12-20
Revised by: Annette Roy, QA Coordinator

Product Name: CytoSep™	**Catalogue No.: M965-40**
This document replaces any previous version	

1. Product Description:

- **Mega Cytology Funnel & Cap:** Disposable, non-sterile funnel for sample volumes of up to 6 ml

2. Packaging:

- **Case:** 12 packages of 40 units / 480 units per case

3. Product Specifications:

Compatible with the Shandon Cytospin® Cytocentrifuge:

- Material:
  - Funnel: Polypropylene
  - Base: Polypropylene
  - Cap: Polyethylene

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **LATEX:** Material is Latex Free
- **CE:** Product is CE marked

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**

1. On exterior case label
2. On label inserted inside the master case
3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination

8. Recommended Use:

- For use in ***Shandon Cytopspin® Cytocentrifuge*** automated equipment, follow the equipment manufacturer's instructions for use

OWNERSHIP OF MATERIALS*:*

Materials and information contained in this document are ***Simport Scientific****'s copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*