VWR We Enable Science

TECHNICAL DATA SHEET

Tap PE, vented, for use with plastic jerricans with external thread (DIN 42).

MATERIAL: Capsule: PE (Polyethylene)
Body: PE (Polyehtylene)


COLOR: Red/Red


Weight: 15.5 ± 1g

Wad: Bore seal

Prod.No.: 29048.098







Remark : Packaging materials and specifications shown may be modified without preceding notice
Version dd. January 27th, 2023

tel/tél +32(0)16 385 011
fax +32(0)16 385 385
info@be.vwr.com
www.vwr.com

Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
RPR/RPM Leuven
BTW/TVA BE 0403.593.343