Nasco SAMPLING
901 Janesville Avenue, Fort Atkinson, WI 53538-0901

Whirl-Pak® Sample Bag
Material Specification Sheet
B01418WA

Size: 4.5" x 9" (11.5 x 23 cm)

Capacity: 18 oz, (532 ml)

Film: Blend of linear low and low-density polyethylene containing black pigmented resins, 2.5 mil (.064 mm) thick, FDA approved for food contact.

Write-On Strip: Horizontal band of thermoset white ink printed on film, compliant with CONEG resolution for heavy metals and European Union Directive for Hazardous Materials.

Closing Tab: Two round aluminum wires covered with plastic tape by a patented process that also covers the wire ends, reduces bag puncture, forms leak proof seal when closed.

Pull Tabs: Made of Tyvek®.

Sterilization: Bags are sterilized with ethylene oxide gas; documentation of sterility is available upon request.

ISO Certification: Manufactured under a quality management system certified to **ISO9001 & ISO13485**.

Packaging: 500 bags/box in a liner bag.

Conditions of Use: Not to be used at temperatures above 180˚F (82˚C).