Dimethylformamide
BAKER ANALYZED® A.C.S. Reagent

Material No.: 9221-03
:
Revision No.: 0

Product Specification

Meets ACS Reagent Chemical Requirements,

Test	Specification
Assay ($HCON(CH_3)_2$) (by GC, corrected for water)	≥ 99.8 %
Appearance	Passes Test
Color (APHA)	≤ 15
Residue after Evaporation	≤ 0.0050 %
Titrable Acid (meq/g)	≤ 0.0005
Titrable Base (meq/g)	≤ 0.003
Water (H_2O)(by Karl Fischer titrn)	≤ 0.1 %
Trace Impurities - Iron (Fe)	≤ 0.1 ppm

For Laboratory,Research,or Manufacturing Use
Meets Reagent Specifications for testing USP/NF monographs

Packaging Site: Paris Mfg Ctr & DC

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700
Avantor Performance Materials, LLC
100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700