
A part of: ThermoFisher SCIENTIFIC

National Scientific
197 Cardiff Valley Rd
Rockwood, TN 37854

Technical Data Sheet

Product Number: 20-SCST
Description: CAPS, PTFE LINED, 20-400, 100/PK

Component	*Specification*
Physical Properties	
Closure Color	White
Description	20-400 Screw Cap
Color	White
Material	Polyethylene Foam, PTFE, Polypropylene
Septum	PTFE/PE Foam Liner
Closure Material	Polypropylene
Quantity	100
thickness (Metric)	1.0 mm
Type	Screw Cap with Septum
Closure Size	20 mm

	01/13/2025	
Revision Number	0	
Approved By	Patrick Giltnane	