888.267.4436
techsupport@origene.com
www.origene.com

Name:Rabbit polyclonal anti-DUS2L antibody
Product Data Sheet - ANTIBODY

Catalog: TA315550

Components:
- Rabbit polyclonal anti-DUS2L antibody (TA315550)

Amount: 100ul

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human DUS2L.

Host: Rabbit

Isotype: IgG

Species Reactivity: Human, Mouse

Guaranteed Applications: WB, IHC

Suggested Dilutions: WB: 1:500~1:3000, IHC: 1:50~1:100, ELISA: 1:20000

Concentration: 1mg/ml

Buffer: Phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Storage Condition: Shipped at -20C. Upon delivery store at -20C. Dilute in PBS (pH7.3) if necessary.
Stable for 12 months from date of receipt.
Avoid repeated freeze-thaws.

Target

Target Name: Homo sapiens dihydrouridine synthase 2 (DUS2), transcript variant 1


Alternative Name: DUS2L; SMM1; URLC8

Database Link: NP_060273


Function:

This product is to be used for laboratory only. Not for diagnostic or therapeutic use.
© 2014 OriGene Technologies, Inc. 9620 Medical Center Dr., Suite 200, Rockville, MD 20850

Validation Data



Western blot analysis of extracts from RAW264.7 cells, using DUS2L antibody.The lane on the right is treated with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human skeletal muscle using DUS2L antibody.The picture on the right is treated with the synthesized peptide.

* More validation images may be available on our website: http://www.origene.com/antibody/TA315550.aspx

This product is to be used for laboratory only. Not for diagnostic or therapeutic use.

© 2014 OriGene Technologies, Inc. 9620 Medical Center Dr., Suite 200, Rockville, MD 20850