interscience

BagClip®

Closing clip

Ref 231 010 (80/100 mL) - **Ref 231 040** (400 mL) - **Ref 231 350** (3500 mL)



Manufactured under:

MADE IN FRANCE

BagClip® allows a temporary closing in 3 points of any type of bag.

TECHNICAL SPECIFICATIONS

- Gas and liquid leakproof
- No contact with the sample
- Available for 80/100/400/2000/3500 mL blender bags
- Designed and made in France

BagClip® 80/100

Ref: 231 010

- Available for 80 and 100 mL bags
- Length (white clip): 130 mm
- Length (green tube): 160 mm
- Box (w x d x h): 21.5 x 15 x 5.5 cm, gross weight: 550 g
- Quantity per pack: 50

BagClip® 400

Ref: 231 040

- Available for 400 mL bags
- Length (white clip): 195 mm
- Length (green tube): 225 mm
- Box (w x d x h): 27 x 19 x 6.5 cm, gross weight: 750 g
- Quantity per pack: 50

BagClip® 3500

Ref: 231 350

- Available for 2000 and 3500 mL bags
- Length (white clip): 395 mm
- Length (green tube): 450 mm
- Box (w x d x h): 51 x 10 x 12 cm, gross weight: 1.6 kg
- Quantity per pack: 50

BagClip®_FT_0321 Pictures and information are not contractually binding. BagClip® is trademark. INTERSCIENCE reserves the right to change or improve the specifications of its products without notice. Please visit www.interscience.com for current updates and additional information. RCS 950 356 220 Versailles. INTERSCIENCE SARL F78860

interscience in the world

PARIS
Phone: +33 1 34 62 62 61
info@interscience.com

FRANKFURT
Phone: +49 611 7238 7770
sales.germany@interscience.com

BOSTON
Phone: +1 781 937 0007
sales.usa@interscience.com

SHANGHAI
电话: +86 178 2123 6642
sales.china@interscience.com

SINGAPORE
Phone: +65 6977 7232
sales.asia@interscience.com

TOKYO
Phone: +81 3 6712 9715
sales.japan@interscience.com

CAUTION

- Do not put in an autoclave
- Do not put in a washing machine
- Can be disinfected with ethanol 70 %, hydrogen peroxyde 6 % or bleach 10 % max

BagClip®_FT_0321 Pictures and information are not contractually binding. BagClip® is trademark. INTERSCIENCE reserves the right to change or improve the specifications of its products without notice. Please visit www.interscience.com for current updates and additional information. RCS 950 356 220 Versailles. INTERSCIENCE SARL F78860

interscience in the world

PARIS
Phone: +33 1 34 62 62 61
info@interscience.com

FRANKFURT
Phone: +49 611 7238 7770
sales.germany@interscience.com

BOSTON
Phone: +1 781 937 0007
sales.usa@interscienceslab.com

SHANGHAI
电话: +86 178 2123 6642
sales.china@interscience.com

SINGAPORE
Phone: +65 6977 7232
sales.asia@interscience.com

TOKYO
Phone: +81 3 6712 9715
sales.japan@interscience.com