Toluene

ChromAR HPLC

Material No: 4483-25
Effective from: 2013/07/11

Product Specification

Test	Specification
Assay (GC)	>= 99.8 %
Water (KF)	<= 0.01 %
Titrable Acid (meq/g)	<= 0.0005
Residue after Evaporation	<= 0.0005 %
Ultraviolet Absorbance (1.00-cm cell vs water)	
290 nm	<= 0.3
300 nm	<= 0.1
310 nm	<= 0.05
330 nm	<= 0.01
350 nm	<= 0.004

For Laboratory, Research or Manufacturing use.
Suitable for Liquid Chromatography and UV-Spectrophotometry

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000