Technical sheet

Brilliant Cresyl Blue Solution

Code product: **310750**

Specifications:

• **Appearance:**	Clear dark blue solution
• **Absorption Maximum in water:**	620 nm ≤ λ ≤ 635 nm
• **Density:**	0.99 - 1.10

Applications:

Haematology: used for examination of reticulocytes on blood smears. Brilliant Cresyl Blue stains reticular filaments deep blue.

Updated: 25/04/2013

RAL Diagnostics - Site Montesquieu - 33651 MARTILLAC Cedex – France - Tel : 33+557 960 404 - Fax : 33+557 960 405 - www.ral-diagnostics.fr