

Product Specifications

Product **444160050**
Ethanol,100%,technical,denaturated with 5% methanol

General Product Data

Version	2
CAS No	64-17-5
Molecular weight	46.06
Molecular formula	C2 H6 O
Linear formula	C2H5OH
Flash point (°C)	12

Product Specifications

Appearance (Color)	Clear colorless
Appearance (Form)	Liquid
Refractive index	1.3590 to 1.3610 (20°C, 589 nm)
Density	(20°C) 0.788 to 0.790 kg/l
GC	>=99.8 % (ethanol) =<5 % (methanol)

Product Specifications are subject to amendment and may change over time.

Acros Organics
ENA23, zone 1, nr 1350, Janssen Pharmaceuticalaan 3a, B-2440 Geel, Belgium
Tel +32 14/57.52.11 - Fax +32 14/59.34.34 Internet: http://www.acros.com
1 Reagent Lane, Fair Lawn, NJ 07410,USA
Fax 201-796-1329