www.abnova.com

9F, No. 108, Jhouzih St.,Taipei, Taiwan
Tel: + 886-2-8751-1888 Fax: + 886-2-6602-1218
E-mail: sales@abnova.com

Datasheet

HSGT1 monoclonal antibody (M01), clone 1A11

Catalog Number: H00011319-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full length recombinant HSGT1.

Clone Name: 1A11

Immunogen: HSGT1 (AAH00721, 1 a.a. ~ 644 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:
MEETMKLATMEDTVEYCLFLIPDESRDSDKHKEILQKYI
ERIITRFAPMLVPYIWQNQPFNLKYKPGKGGVPAHMF
GVTKFGDNIEDEWFIVYVIKQITKEFPELVARIEDNDGE
FLLIEAADFLPKWLDPENSTNRVFFCHGELCIIPAPRKS
GAESWLPTTPPTIPQALNIITAHSEKILASESIRAAVNRR
IRGYPEKIQASLHRAHCFLPAGIVAVLKQRPRLVAAAV
QAFYLRDPIDLRACRVFKTFLPETRIMTSVTFTKCLYAQ
LVQQRFVPDRRSGYRLPPPSDPQYRAHELGMKLAHG
FEILCSKCSPHFSDCKKSLVTASPLWASFLESLKKNDY
FKGLIEGSAQYRERLEMAENYFQLSVDWPESSLAMSP
GEEILTLLQTIPFDIEDLKKEAANLPPEDDDQWLDLSPD
QLDQLLQEAVGKKESESVSKEEKEQNYDLTEVSESMK
AFISKVSTHKGAELPREPSEAPITFDADSFLNYFDKILG
PRPNESDSDDLDDEDFECLDSDDDLDFETHEPGEEAS
LKGTLDNLKSYMAQMDQELAHTCISKSFTTRNQVEPV
SQTTDNNSDEEDSGTGESVMAPVDVDLNLVSNILESY
SSQAGLAGPASNLLQSMGVQLPDNTDHRPTSKPTKN

Host: Mouse

Reactivity: Human

Applications: ELISA, IF, S-ELISA, WB-Ce, WB-Re
(See our web site product page for detailed applications information)

Protocols: See our web site at http://www.abnova.com/support/protocols.asp or product page for detailed protocols

Isotype: IgG1 lambda

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 11319

Gene Symbol: ECD

Gene Alias: GCR2, HSGT1