Specification Sheet

Product Name	Quinoline Yellow 2SF analytical standard
Product Number	01354
Product Brand	SIAL
CAS Number	83-08-9
Molecular Weight	273.29

TEST	SPECIFICATION
APPEARANCE (COLOR)	Yellow to Dark Yellow and Yellow-Orange
APPEARANCE (FORM)	Powder or Crystals
PURITY (HPLC AREA %)	≥ 97.0 %
IDENTIFICATION (LC-MS)	CORRESPONDS TO REQUIREMENTS
INFRARED SPECTRUM	CONFORMS TO STRUCTURE

TEST	SPECIFICATION
EXPIRY PERIOD	36 MONTHS