
CATALOG NUMBER	SEE BELOW	
WORKING RANGE (µL)	50 – 250µL	
MATERIAL	PS	
	MM	INCH
A, BASE LENGTH	127,7	5,03
B, BASE WIDTH	85,7	3,37
C, TOP LENGTH	123,3	4,85
D, TOP WIDTH	81,2	3,20
E, FLANGE LONG SIDE	2,7	0,11
F, OVERALL HEIGHT	14,5	0,57
G1, DIA TOP	Ø6,8	Ø0,27
G2, DIA BOTTOM	Ø6,4	Ø0,25
H, A - I LOCATION (X)	14,4	0,57
J, A - I LOCATION (Y)	11,4	0,45
K, WELL TO WELL SPACING	9	0,35
L, WELL DEPTH	10,3	0,41
M, BOTTOM OF WELL DISTANCE	4,3	0,17
N, FLANGE SHORT SIDE	2,7	0,11

CATALOG NO	SURFACE	COLOR	STERILE	WITH LID*	UNITS/ PK/CASE
249662	NON-TREATED	CLEAR	+	-	1/50

NOTES:
1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATION PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

PLATE

SECTION A-A



ThermoFisher SCIENTIFIC

Roskilde Site · Kamstrupvej 90 · P.O.Box 280
DK-4000 Roskilde · Denmark
Tel +45 4631 2000 · Fax +45 4631 2042
E-mail: info.nunc@thermofisher.com

www.nuncbrand.com

SUBJECT:	DRAWN BY/DATE:	CHANGE ORDER NO.::	VALID FROM:	APPROVED BY/DATE:	FORMA
V96 MICROWELL PLATE, PS					
	ITEM NUMBER:	MATERIAL:	DRAWING NO.:	VERSION:	PAGE TYPE: SCALE: METRI

CONFIDENTIAL
THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETORY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC