
Technical Data Sheet





VWR® Microscope Slides

Description:	Clear imprinted with VWR LOGO ; Clear with a “●” at each of the two lower corners; pre-cleaned, clear; no marks, prints or contamination, etc. Color: Standard colors Color Frosted Area: 19.5 x 25 (±0.5)mm;
Purpose:	Used for histological pathology examination
Materials:	Standard Glass; thinckness 1.0 (-0.05,+0.1)mm
Features:	The smooth edges provide the safety of the operation
	Color Frosted
	Suitable for marking with conventional labeling system, pencil, mark pens and thermal transfer printers
Temperature	Normal laboratory temperature and humidity
Uses:	Recommended for routine H&E stains,frozen sections
Shelf life :	60 Months

VWR® Microscope Slides

EU Part No.	Dimension				Frosted or Painted	Paint Color	Corner Angle	Edges	Packaging	QTY/CASE
	A (mm)	B (mm)	C (mm)	D (mm)						
			-0.05/+0.1	(+/-0.5)						
631-1169	75.7~76.0	25.0~25.3	1	N/A	Painted	BLUE	45°	Smooth	50pcs / Box, 50 Boxes / Case	2500
631-1172	75.7~76.0	25.0~25.3	1	N/A	Painted	GREEN	45°	Smooth	50pcs / Box, 50 Boxes / Case	2500
631-1173	75.7~76.0	25.0~25.3	1	N/A	Painted	ORANGE	45°	Smooth	50pcs / Box, 50 Boxes / Case	2500
631-1175	75.7~76.0	25.0~25.3	1	N/A	Painted	PINK	45°	Smooth	50pcs / Box, 50 Boxes / Case	2500
631-1176	75.7~76.0	25.0~25.3	1	N/A	Painted	YELLOW	45°	Smooth	50pcs / Box, 50 Boxes / Case	2500
631-1183	75.7~76.0	25.0~25.3	1	N/A	Painted	WHITE	45°	Smooth	50pcs / Box, 50 Boxes / Case	2500