Industrial bucket with metal handle



Industrial bucket for transporting, decanting and mixing. Safe emptying with wide spout. Handle made of metal.

- LLDPE, white
- With litre scale, graduation 1000 ml
- With spout

Industrial bucket with metal handle				
Capacity l	Dia. mm	Height mm	Item no.	
9	280	250	2308-0900	N3
12	300	280	2308-1200	N3
17	340	310	2308-1700	N3

Laboratory bucket



Laboratory bucket for transporting, decanting and mixing. Safe emptying with wide spout. Handle made of metal.

- PE, white
- With contents scale, graduation 1000 ml
- With spout

Laboratory bucket				
Capacity l	Dia. mm	Height mm	Item no.	
10	290	280	2306-0010	N3
15	320	340	2306-0015	N3