

Product description

Order number	95.64.997
Product description	Storage box, slip-on lid, material: cardboard, white, format: 10 x 10, for 100 collection tubes, suitable for tubes up to max. 12 mm Ø, tube height between 36-45 mm, 1 piece(s)/bag

Product characteristics

Closure type	slip-on lid
Suitable for	tubes up to max. 12 mm Ø, tube height between 36-45 mm
Storage capacity	100
Format	10 x 10

Size

Width of product	135 mm
Height of product	45 mm
Length of product	135 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Material & colours

Product material	cardboard
Colour of product	white
Closure material	cardboard
Colour of cap	white

Purity & certification

Product category	no medical device | IVD
Batched	yes

Packaging

Minimum order qty.	1
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	120
Piece(s) / pallet	1440
Depth of case	546 mm
Width of case	419 mm
Height of case	495 mm
Case volume	0.1134 cbm
Weight of product	0.0875 kg
Weight of case	11.5 kg
EAN of inner box	4038917095907
EAN case	4038917007887

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.