


ORIGINAL

Tolerances unless otherwise specified:
Without decimal: ±0.5 mm
1 place decimal: ±0.1 mm
2 place decimal: ±0.05 mm

Undefined radius: R 0.20 mm

28/2-01 PCØ

SUBJECT: Loop 10 micro liter. 251586		SCALE: 2:1(10:1)	COMPLETE BY/DATE: TBS/24.01.00
		VERSION: 1.1	APPROVED BY/DATE: 2/24.02.00
		MATERIAL: PS	VALID FROM: 24.02.00
		DRAWING: 2158	PAGE 1 OF 1

Nunc A/S
Kornstrupvej 90
P.O. Box 280
Dk-4000 Roskilde
Denmark
Tel 45 46 35 90 65
Fax 45 46 32 08 99
e-mail: infociety@nunc.dk
A Subsidiary of Sybron International

28/02/01 14:10 NUNC 4000 ROSKILDE DK → NUNC INC 1 NR.677 S006/007