A part of: ThermoFisher SCIENTIFIC

National Scientific
197 Cardiff Valley Rd
Rockwood, TN 37854

Technical Data Sheet

Product Number: C4011-88
Description: 2ml Target Clear Vial/Seal 12x32mm 500/pk

Component	***Specification***
Vial Physical Properties	
Material	Clear borosilicate glass
Height	32.0mm
Outside Diameter	11.65mm
Inside Opening Diameter	6.25mm
Wall Weight	0.8mm
Usable capacity	1.5ml
Cap	
Cap material	Aluminum
Type	Crimp Seal
Cap Color	Natural Silver
Opening Diameter	5.3mm
Septa Material	PTFE/ Natural Rubber

Date	11/3/23	
Revision Number	0	
Approved By	Patrick Giltnane	