interscience

# BagLight® PolySilk®

## Blender bag without filter, transparent

**Ref. 131 025** (100 mL) – **Ref. 132 025** (400 mL) – **Ref. 132 050** (400 mL) – **Ref. 132 200** (2000 mL)
**Ref. 133 025** (3500 mL)



Standards:

ISO 7218 | ISO 6887 | FDA BAM Bacteriological Analytical Manual

Manufactured under:

MADE IN FRANCE

**BagLight® PolySilk®** is a blender bag without filter. It is ideal for sample preparation.

### TECHNICAL SPECIFICATIONS

- Non-filter bag
- **PolySilk®**, polyolefin complex
- Flexible and transparent
- Resistant to freezing and high temperatures (from -40°C to 80°C / -40°F to 176°F)
- No contact between the sample and the blender during homogenization of the sample
- Compatible with any lab blender
- Gamma ray treated: Gamma 5 to 12 kGy, with certificate
- Shelf-life: 10 years
- Storage: in a temperate and dry place, away from humidity and light
- Available in 100/400/2000/3500 mL

- In compliance with ISO 7218, ISO 6887 and FDA BAM (Bacteriological Analytical Manual)
- Approved for food contact - Regulation (EC) N° 1935/2004
- Designed and made in France

## BagLight® PolySilk® 100

**Ref. 131 025**

- Max blending volume: 100 mL
- Optimal blending volume: 5 - 50 mL
- Bag dimensions: 110 x 200 mm
- Box dimensions: 31 x 22 x 7 cm, weight: 1.80 kg
- Pack of 25
- Box of 500

## BagLight® PolySilk® 400

**Ref. 132 025**

- Max blending volume: 400 mL
- Optimal blending volume: 50 - 300 mL
- Bag dimensions: 175 x 300 mm
- Box dimensions: 35 x 21 x 11 cm, weight: 4.10 kg
- Pack of 25
- Box of 500



**BagLight® PolySilk** _FT_0322 Pictures and information are not contractually binding. BagLight® is a trademark. INTERSCIENCE reserves the right to change or improve the specifications of its products without notice. Please visit www.interscience.com for current updates and additional information. RCS 950 356 220 Versailles. INTERSCIENCE SARL F78860

interscience in the world

PARIS
Phone: +33 1 34 62 62 61
info@interscience.com

FRANKFURT
Phone: +49 611 7238 7770
sales.germany@interscience.com

BOSTON
Phone: +1 781 937 0007
sales.usa@interscience.com

SHANGHAI
电话: +86 178 2123 6642
sales.china@interscience.com

SINGAPORE
Phone: +65 6977 7232
sales.asia@interscience.com

TOKYO
Phone: +81 3 6712 9715
sales.japan@interscience.com

## BagLight® PolySilk® 400

**Ref. 132 050**

- Max blending volume: 400 mL
- Optimal blending volume: 50 - 300 mL
- Bag dimensions: 175 x 300 mm
- Box dimensions: 35 x 21 x 11 cm, weight: 4.10 kg
- Pack of 50
- Box of 500

## BagLight® PolySilk® 2000

**Ref. 132 200**

- Max blending volume: 2000 mL
- Optimal blending volume: 400 - 1500 mL
- Bag dimensions: 250 x 380 mm
- Box dimensions: 48 x 29 x 14 cm, weight: 10.2 kg
- Pack of 25
- Box of 500

## BagLight® PolySilk® 3500

**Ref. 133 025**

- Max blending volume: 3750 mL
- Optimal blending volume: 400 - 3750 mL
- Bag dimensions: 380 x 560 mm
- Box dimensions: 41 x 34 x 12 cm, weight: 7.70 kg
- Pack of 25
- Box of 250



**BagLIght® PolySilk** _FT_0322 Pictures and information are not contractually binding. BagLight® is a trademark. INTERSCIENCE reserves the right to change or improve the specifications of its products without notice. Please visit www.interscience.com for current updates and additional information. RCS 950 356 220 Versailles. INTERSCIENCE SARL F78860

interscience in the world

PARIS
Phone: +33 1 34 62 62 61
info@interscience.com

FRANKFURT
Phone: +49 611 7238 7770
sales.germany@interscience.com

BOSTON
Phone: +1 781 937 0007
sales.usa@intersciencelab.com

SHANGHAI
电话: +86 178 2123 6642
sales.china@interscience.com

SINGAPORE
Phone: +65 6977 7232
sales.asia@interscience.com

TOKYO
Phone: +81 3 6712 9715
sales.japan@interscience.com