LABORATORY WARE



BEAKERS MOLYBDENE STAINLESS STEEL

Ref.	Diam.	Height	Capacity
Without spout			
2010	50 mm	60 mm	0,100 L
2025	65 mm	80 mm	0,250 L
2050	80 mm	110 mm	0.500 L
2100	100 mm	130 mm	1 L
2200	120 mm	180 mm	2 L
2300	150 mm	180 mm	3 L
2500	180 mm	210 mm	5 L

For option ordering please add after the reference **B** for spout.
A for handle, **AB** for handle and spout.

HIGH EVAPORING DISHES

Ref.	Diam.	Height	Capacity
1070	70 mm	35 mm	125 ml
1085	85 mm	45 mm	250 ml
1110	110 mm	55 mm	500 ml
1140	140 mm	70 mm	1 L
1160	160 mm	105 mm	2 L
1200	200 mm	160 mm	5 L

CEVAPORING DISHES

Ref.	Diam.	Height	Capacity
406	60 mm	20 mm	50 ml
408	80 mm	20 mm	100 ml

FUNNELS WITH HANDLE

Ref.	Diam.	Capacity
4215	150 mm	0,5 litreL
4220	195 mm	1 L
4250	250 mm	3 L

SPATULAS BY 5

Réf.	Model	Length
1869	Stainless steel spatula	195 mm
1975	Double spatula	230 mm
1822	Shaft spatula	130 mm
1826	Micro-spatula	130 mm
1975LG	Double spatula	500 mm

CHEMICAL SCOOP

Ref.	Diam.	Height	Capacity
915	75 mm	150 mm	0,5 L
920	95 mm	200 mm	1 L
923	115 mm	233 mm	2 L