
PROSAT® Theta™ Wipes

These PROSAT® Theta™ wipes are manufactured from a lightweight hydroentangled cellulose/polyester and saturated with a solution of 70% USP grade IPA and 30% deionized water. A cost effective wipe ideal for transfer disinfection and general cleaning, PROSAT Theta wipes offer a soft and durable wipe providing consistant solvent application.

Presaturated wipes ensure consistent saturation of each wipe independent of operator. Presaturated wipes can increase solvent control and accountability as well as reduce VOC emissions. The wipes are provided in convenient and easy to use peel and reseal pouches. The resealable pouches preserve solvent saturation and cleanliness.



Features	Benefits
Nonwoven material	• Low in particles and fibers making an excellent general purpose wipes
Presaturated wipes	• To reduce solvent usage and VOC emissions • Ensures consistent saturation of each wipe independent of the operator
Resealable pouch	• Wipes are provided in a resealable pouch for ease of use and convenience
Validated sterile via gamma irradiation to a 10^{-6} Sterility Assurance Level	• Suitable for use in Grade A and B cleanrooms

Part No.	Description	Packaging
PSC20005	PROSAT® *Sterile*™ Theta Presaturated with 70% IPA and 30% DI water, 9" x 11" (23x28cm)	50/pouch; 28 pouches/case
PSC20006	PROSAT® Theta Presaturated with 70% IPA and 30% DI water, 9" x 11" (23x28cm)	50/pouch; 28 pouches/case
PSCP0001	PROSAT® *Sterile*™ Theta Presaturated with 70% IPA and 30% DI water, 6" x 9" (15x23cm)	50/pouch; 40 pouches/case

Copyright © 2022 Contec, Inc. All rights reserved.

Contec, Inc
Spartanburg, SC
United States
Tel: +1 (864) 503-8333


Contec
Vannes Cedex
France
Tel: +33(0) 2 97 43 76 98

Contec Cleanroom Technology (Suzhou) Co. Ltd
Suzhou
China
Tel: +86-512-6274 4050

Product Information

Material	54% Cellulose/ 46% polyester
Construction	Hydroentangled staple fibers
Packaging Materials	Pouch (PCH), low density polyethylene (LDPE) and polyester (PET) Flow-Wrap Outer Bag (FOB), low density polyethylene (LDPE) and polyester (PET) Outer Bag (OB1, OB2), low density polyethylene (LDPE) Case (CS), corrugated fiberboard (PAP)
Solution of saturation	70% USP Grade isopropyl alcohol/ 30% deionized water
Environment	ISO 5-8 Grade A/B for sterile, C/D for nonsterile





Technical Data

Attribute (units)	Typical Value	Test Method
Basis weight; nominal (g/m^2)	53	Contec Method
Non-volatile residue, NVR		IEST-RP-CC004.3, Sec. 7.1.2
In deionized water; (g/m^2)	0.037	
In isopropyl alcohol; (g/m^2)	ND	
Specific ions		IEST-RP-CC004.3, Sec. 7.2.2
Sodium; (ppm)	11.9	
Chloride; (ppm)	18.2	
Particles, readily releasable		IEST-RP-CC004.2, Sec. 5.1
$P \geq 0.5\mu m$; ($x10^6/m^2$)	18.6	
Fibers > 100µm; ($x10^3/m^2$)	21.7	

Recycle Symbol Key

PET	1
HDPE	2
LDPE	4
PP	5
PAP	20

Packaging

	EA/PCH	PCH/OB1	OB1/CS	EA/CS
PSC20006	50	4	7	1400

	EA/PCH	PCH/FOB	FOB/OB2	OB2/CS	EA/CS
PSC20005	50	1	7	4	1400
PSCP0001	50	1	8	5	2000

EA = Each, FOB = Flow-Wrap Outer Bag, OB1 = Outer Bag 1, OB2 = Outer Bag 2, PCH = Pouch, CS = Case

Notes

a) The data shown are typical values and should not be used as product specifications.
b) Valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
c) Current and/or comparison data may be available. Please contact a Contec sales representative for details.

Test Methods:

1) CTM Contec Test Method
2) IEST-RP-CC004.3 Evaluating Wiping Materials Used in Cleanroom and Other Controlled Environments, Institute of environmental Sciences and Technology, Rolling Meadows IL.

Contec, Inc
Spartanburg, SC
United States
Tel: +1 (864) 503-8333

Contec
Vannes Cedex
France
Tel: +33(0) 2 97 43 76 98

Contec Cleanroom Technology (Suzhou) Co. Ltd
Suzhou
China
Tel: +86-512-6274 4050

www.contecinc.com