Product Specifications Sheet

VWR WOMEN´S HACCP COAT WITH STAND-UP COLLAR

Coat women´s style with stand-up collar made from 65% Polyester/35% Cotton, 2 inside side-pockets, 1 inside breast-pocket on the left, long sleeves with 2-step adjustable cuffs with 2 V2A-steel press studs , 2 side slits, concealed press stud band with 5 V2A-steel press studs, fitted, with back belt

Intended use

Food handling
Food manufacturing areas

Specifications

Material	65% Polyester/35% Cotton, DIN 10524
Weight (approx.)	215 g/m²
Color	White
Fibre characteristic	OEKO-TEX Standard 100
Washing instruction	 suitable for industrial laundering



VWR International bvba/sprl Researchpark Haasrode 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium
Tel.: 016 385 011 • Fax: 016 385 385 • E-mail: info@be.vwr.com

Version 1:Date of issue: **May, 2012**

Product Specifications Sheet

VWR WOMEN´S HACCP COAT WITH STAND-UP COLLAR

Sizes				Measurements in cm				Catalogue number
D	IT	UK	FR	½ Chest size (1)	Crown width (2)	Total (3)	Arm length (incl. cuffs) (4)	
34	40	8	36	47,5	37,9	95,5	57	113-1804
36	42	10	38	49,5	39,3	97	58	113-1805
38	44	12	40	51,5	40,7	98,5	59	113-1806
40	46	14	42	53,5	42,1	100	60	113-1807
42	48	16	44	55,5	43,5	101,5	61	113-1808
44	50	18	46	57,5	44,9	103	62	113-1809
46	52	20	48	59,5	46,3	104,5	63	113-1810
48	54	22	50	61,5	47,7	106	64	113-1811
50	56	24	52	63,5	49,1	107,5	65	113-1812
52	58	26	54	65,5	50,5	109	66	113-1813



Product compliance:

Council Directive 89/686/EEC relating to Personal protective equipment, CAT I against low risks, DIN 10524 Risk class 1-3
Certificate of conformity and DIN 10524 available

Disclaimer: *VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.*

VWR International bvba/sprl Researchpark Haasrode 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium
Tel.: 016 385 011 • Fax: 016 385 385 • E-mail: info@be.vwr.com

Version 1:Date of issue: May, 2012