BAKERBOND spe, NARC-2 Carbon

BAKERBOND spe

Disposable Extraction Columns

Material No: 7225-05
Effective from: 2015/06/07

# Product Specification

| Test | Specification |
|---|---|
| The following results are based on the bulk bonded phase: | |
| Carbon (C)(%) | For Information Only |
| Sulfur (S)(ug/ml) | For Information Only |
| Surface Coverage (ug/m2), Based on %C | For Information Only |
| Surface Coverage (ug/m2), Based on %S | For Information Only |
| The following results are based on the silica gel raw material: | |
| Particle Size Distribution: APD, microns | 47 - 61 |
| Specific Surface Area, m2/g | For Information Only |
| The following results are based on the solid phase extraction columns: | |
| Extraction Residue | <= 0.2 % |

For Laboratory, Research or Manufacturing use.
6mL Solid Phase Extraction Columns
250 mg per Column

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000