BAKERBOND spe, SDB

BAKERBOND spe

Disposable Extraction Columns; Glass

Material No: 7609-02
Effective from: 2023/11/08

# Product Specification

| Test | Specification |
|---|---|
| Styrene Divinylbenzene (SDB-1): | |
| Breakthrough Antrazine* | >= 1 L |
| Breakthrough Desisopropylatrazine* | >= 1 L |
| Specific Surface Area | >= 950 m2/g |
| Porosity, Median Pore Diameter | 50 - 70 A |
| Particle Size - 50% Volume | 65 - 95 um |
| Capacity (mg caffeine / g sorbent) | >= 250 |
| Recovery Atrazine* | >= 97 % |
| Recovery Desisopropylatrazine* | >= 75 % |

For Laboratory, Research or Manufacturing use.
3mL Solid Phase Extraction Columns
200 mg per Column

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000