Seit 1927
über 90 Jahre

Since 1927
more than 90 Years

Amarell
GmbH & Co. KG
Präzisions-Thermometer / Aräometer

Lindenstraße 43
D-97892 Kreuzwertheim

Tel. +49 (0) 9342 9283-0
Fax: +49 (0) 9342 39860

Internet: www.amarell.de
E-Mail: thermometer@amarell.de

Technical Data Sheet

Article number:	L25992-FL
Type:	Precision thermometer, enclosed scale, capillary colourless
Measurement range:	0+50°C
Subdivision:	0,1°C
Length:	420mm
Diameter:	8,5-9,0mm
Filling:	Red special liquid
Immersion:	Total
Finish:	Richter-top finish with button

March‘25