PRODUCT DATA SHEET

Injekt® (2 ml, luer, non-sterile) (9202528)

2-piece single-use syringe

Attribute Name	Attribute Value
Product type	Syringe
Application	Standard
Construction of syringe	2-piece
Style of syringe cone	Luer
Position of syringe cone	centric
Nominal capacity	2.0 ml
Total graduated capacity	3.0 ml
Max. usable capacity	3.30 ml
Scale value	0.1000 ml
Light-protected	No
Color of light protection	-
Type of packaging	Multipack
Packaging contains PVC	No
Packaging contains Latex	No

B. Braun Melsungen AG
Carl-Braun-Straße 1
34212 Melsungen | Germany
Tel. +49 5661 710
www.bbraun.com

Attribute Name	Attribute Value
Packaging contains DEHP	No
Product (unpack.) contains PVC	No
Product (unpack.) cont. Latex	No
Product (unpack.) cont. DEHP	No
Product contains PVC	No
Product contains Latex	No
Product contains DEHP	No
Product contains Nickel	No
Product contains BPA	No
Sterility	non-sterile
Type of Sterilization	not sterilized ex works
Type of Sterilizat. (allowed)	Unsteril
Disposable / Reusable	disposable