N,N-Dimethylformamide

ChromAR HPLC

Material No: 5362-25
Effective from: 2013/07/09

Product Specification

Test	Specification
Assay (GC)	>= 99.8 %
Water (KF)	<= 0.05 %
Titrable Acid (meq/g)	<= 0.0005
Residue after Evaporation	<= 0.0005 %
Ultraviolet Absorbance (1.00-cm cell vs water)	
275 nm	<= 0.2
280 nm	<= 0.1
290 nm	<= 0.06
300 nm	<= 0.03
320 nm	<= 0.01

For Laboratory, Research or Manufacturing use.
Suitable for Liquid Chromatography and UV-Spectrophotometry

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000