


SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETORY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	APP'D		REV.
PB	11/29/07	PSS 00-56X-0020-1	1107
CK'D			

DESCRIPTION:

aPES FILTER UNITS


CAT. NO. . . . . :	565-0020		568-0020		566-0020	
UPPER CAP. (mL). . . . :	150		250		500	
RECEIVER CAP. (mL). :	150		250		500	
PORE SIZE (um). . . . . :	0.2		0.2		0.2	
MEMBRANE COLOR. . . . . :	WHITE		WHITE		WHITE	
MEMBRANE DIA (mm). . :	50		50		75	
	INCHES	MM	INCHES	MM	INCHES	MM
A (UPPER I.D.). . . . . :	3.15	80.0	3.15	80.0	3.80	96.5
B HEIGHT. . . . . :	6.15	156.2	7.9	200	9.8	249
C DIA (BASE). . . . . :	3.50	88.9	3.50	88.9	3.86	98.0
D OVERALL. . . . . :	3.7	94	3.7	94	4.0	102
NOM. WT. (gms). . . . . :	78		95		155	

MATERIALS	
UPPER CONTAINER	PS
COVER	PS
COLLAR	BLUE PS
RECEIVER	PS
STORAGE CLOSURE	WHITE HDPE
ADAPTOR	WHITE LDPE
GASKET	FOAMED LDPE
MEMBRANE	aPES
VENT PLUG	CA

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CAT. NO.	569-0020		567-0020	
PORE SIZE (um)	0.2		0.2	
UPPER CAP. (mL)	500		1000	
RECEIVER CAP. (mL)	500		1000	
MEMBRANE COLOR	WHITE		WHITE	
MEMBRANE DIA (mm)	90		90	
	INCHES	MM	INCHES	MM
A (UPPER I.D.)	3.84	97.5	4.48	113.8
B HEIGHT	10.7	272	12.8	325
C DIA (BASE)	3.86	98.0	4.48	113.8
D OVERALL	4.0	102	4.2	106
NOM. WT. (gms)	170		250	
MATERIALS				
UPPER CONTAINER	PS			
COVER	PS			
COLLAR	BLUE PS			
RECEIVER	PS			
STORAGE CLOSURE	WHITE HDPE			
ADAPTOR	WHITE LDPE			
GASKET	FOAMED LDPE			
MEMBRANE	aPES			
VENT PLUG	CA			



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETORY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.


ThermoFisher SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	BR 5/6/10	PSS 00-056X-0020-2	REV 0510
CK'D	BR				

DESCRIPTION

MEDIA-PLUS aPES FILTER UNITS



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETORY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	APP'D		REV
PB		PSS-564-0020	1107
CK'D	11/29/07		

DESCRIPTION:

FILTER UNIT - 50mL CAPACITY

CAT. NO.....:	564-0020	
UPPER CAP (mL).....:	150	
TUBE CAP (mL).....:	50	
PORE SIZE (µm).....:	0.2	
MEMBRANE COLOR....:	WHITE	
MEMBRANE DIA (mm):	50	
CLOSURE SIZE (mm)....:	33	
	INCH	MM
A (UPPER I.D.).....:	3.16	80.3
B HEIGHT.....:	7.40	188.0
C DIA.....:	4.00	101.6
NOM. WT. (gms).....:	75	

MATERIALS	
COVER	PS
UPPER CONTAINER	PS
ADAPTOR (FILTER)	WHITE LDPE
COLLAR	BLUE HIPS
ADAPTOR TUBE/COLLAR	BLUE HIPS
TUBE	PP
STAND	WHITE HIPS
CLOSURE	HDPE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.

2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.