





REF European Cat no	Product Description	Capacity (ml)	Brim Volume (ml)	Material (Bottle / Cap)	Bottle Shape	Colour (Bottle / Cap)	Neck Size & Finish	Neck ID (mm)	Length (mm)	Width (mm)	Height Without Cap (mm)	Height With Cap (mm)	General Wall Thickness (mm)	Bottle Weight (Grm)	Cap Weight (Grm)	Gross Weight (Grm) ±0.5grm
215-4247	Urine Collection Bottle, square, red cap	150	172	HDPE / PP	Square	Natural / RED	45 / 400	37.42	50.5	50.5	82.24	86.2	0.7mm	16	4.5	20.5
215-4248	Urine Collection Bottle, square, red cap	250	272	HDPE / PP	Square	Natural / RED	45 / 400	37.42	57	57	100.24	102.2	0.7mm	20	4.5	24.5
215-4249	Urine Collection Bottle, square, blue cap	250	272	HDPE / PP	Square	Natural / BLUE	45 / 400	37.42	57	57	100.24	102.2	0.7mm	20	4.5	24.5
215-4250	Urine Collection Bottle, square, red cap	500	550	HDPE / PP	Square	Natural / RED	63 / 410	56.12	76.5	76.5	119.5	122.5	0.7mm	40	10.5	50.5
215-4251	Urine Collection Bottle, square, blue cap	500	550	HDPE / PP	Square	Natural / BLUE	63 / 410	56.12	76.5	76.5	119.5	122.5	0.7mm	40	10.5	50.5
215-4252	Urine Collection Bottle, square, white cap	500	550	HDPE / PP	Square	Natural / WHITE	63 / 410	56.12	76.5	76.5	119.5	122.5	0.7mm	40	10.5	50.5
215-4253	Urine Collection Bottle, rectangular, red cap	1000	1070	HDPE / PP	Rectangular	Natural / RED	63 / 410	56.12	94.9	76.3	178.2	181.18	0.7mm	55	10	65
215-4254	Urine Collection Bottle, rectangular, blue cap	1000	1070	HDPE / PP	Rectangular	Natural / BLUE	63 / 410	56.12	94.9	76.3	178.2	181.18	0.7mm	55	10	65
215-4255	Urine Collection Bottle, rectangular, white cap	1000	1070	HDPE / PP	Rectangular	Natural / WHITE	63 / 410	56.12	94.9	76.3	178.2	181.18	0.7mm	55	10	65
215-4256	Urine Collection Bottle, rectangular, red cap	2000	2120	HDPE / PP	Rectangular	Natural / RED	63 / 410	56.12	108	97	230	233	0.7mm	90	10	100
215-4257	Urine Collection Bottle, rectangular, white cap	2000	2120	HDPE / PP	Rectangular	Natural / WHITE	63 / 410	56.12	108	97	230	233	0.7mm	90	10	100
215-4258	Urine Collection Bottle, square, white cap	2500	2700	HDPE / PP	Square	Natural / WHITE	83 / 2000	71.5	109	109	257	264.5	0.7mm	112	18	130
215-4259	Urine Collection Bottle, amber square, yellow Cap	2700	3100	HDPE / PP	Rectangular	Amber / Yellow	83 / 2000	71.5	125	115	250	253.36	0.7mm	122	18	140
215-4361	Urine Collection Bottle, square with plug & red Cap	500	550	HDPE / PP	Square	Natural / RED	63 / 410	56.12	76.5	76.5	119.5	122.5	0.7mm	40	10.5	50.5
215-4362	Urine Collection Bottle, rectangular with plug & red cap	1000	1070	HDPE / PP	Rectangular	Natural / RED	63 / 410	56.12	94.9	76.3	178.2	181.18	0.7mm	55	10	65
215-4363	Urine Collection Bottle, rectangular with plug & red cap	2000	2120	HDPE / PP	Rectangular	Natural / RED	63 / 410	56.12	108	97	230	233	0.7mm	90	10	100

Material

Bottle : HDPE, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed.
Closure : Polypropylene, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed With Food grade colour.

BSE/TSE Statement: These product is manufactured without the use of Raw materials of animal origin.

Latex Free Statement: No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.



Material Recycle Code

2 HDPE 4 LDPE 5 PP

Description: Urine Collection Bottle

Drawing: Not to Scale

VWR

Release: 2020/08/27

Revision: 0