:\GCA\CSD-635525


TH IS DRAWING CONTAINS CONFIDENTIAL INFORMATION AND IS A COPYRIGHT OF DWK LIFE SCIENCES INC. IT MAY NOT BE DUPLICATED FOR ANY PURPOSE IN WHOLE OR IN PART WITHOUT THE EXPRESS WRITTEN CONSENT OF DWK LIFE SCIENCES INC. ALL RIGHTS RESERVED.©
ELECTRONIC OR PAPER COPIES OF THIS PRINT OUTSIDE OF DWK LIFE SCIENCES SYSTEM ARE UNCONTROLLED AND ARE FOR REFERENCE ONLY
UNLESS OTHERWISE STAMPED.

TOLERANCE: ± 2mm UNLESS NOTED

THIRD ANGLE PROJECTION

PRODUCT COMPONENTS

ITEM	QTY	DESCRIPTION



MATERIAL: TYPE I CLEAR BOROSILICATE GLASS, 33 EXPANSION.

A	TITLE 1000mL FLASK, FILTER, 40/35 JT	DATE DRAWN 05-11-90	PROJ#:	–	REVISED TO DWK STDS	12/09/19 EL
		DRAWN BY K.DAVIDSON	APPROX. WGT.	A	12370	04-03-01 DB
	CAPACITY 1,000 mL	CHECKED BY T.R.J.	AUTHORIZED R.B.	REV	E.C.N. #	DATE
1. DIMENSIONS ARE IN MILLIMETERS.	FINISH ⅀ 40/35 GND INNER JOINT	DURAN WHEATON KIMBLE 1501 North 10th Street | Millville, NJ 08332-2038 | USA		DWG. NO: B- CSD-635525		
----NOTES----	SCALE: FULL UNLESS NOTED. DO NOT SCALE THIS DRAWING.					