MICROSOLV TECHNOLOGY CORPORATION

Leaders in Separation & Purification Technologies

Product Data Sheet

9509S-3DAP-PX

Vial & Cap kit. 100 each 2 mL, screw top, hydrophobic, amber autosampler vials and matching black caps with clear, pre-slit, silicone rubber / PTFE septa.



RSA-ProX™ Brand

- Vials treated with RSA-Pro X™ technology to prevent their already pristine surface from adsorbing stubbornly sticky compounds that have affinity for the glass structure. The Hydrophobic surface is not only non-adsorbing for most compounds but is extremely stable even in the presence of water.
- The hydrolytic stability will last for years compared to Silanized Vials.
- Amber (low actinic glass) is used for light sensitive compounds.
- AQR brand caps, black with inert clear, silicone rubber/PTFE fitted septa.
- 2 mL nominal volume.
- 12 x 32 mm outer dimensions of the vial with a 9-425 mm screw thread top.
- Caps have matching 9-425 mm threads.
- This product is manufactured in an ISO 9001:2015 Certified facility.

9158 Industrial Blvd NE
Leland, NC 28451

p: 1.732.380.8900
f: 1.910.769.9435

info@mtc-usa.com
www.mtc-usa.com