Methanol

Baker HPLC Analyzed

Isocratic Grade

Material No: 8404.9030RC
Effective from: 2015/12/20

Product Specification

Test	Specification
Assay	>= 99.8 %
Color (APHA)	<= 10
Residue after Evaporation (in ppm)	<= 5
Titrable Acid (meq/g)	<= 0.0003
Titrable Base (meq/g)	<= 0.0002
Water (H2O)	<= 0.05 %
Ultraviolet Absorbance (1.00-cm path vs water):	
at 225 nm	<= 0.30
at 240 nm	<= 0.10
at 265 nm	<= 0.01

For Laboratory, Research or Manufacturing use.
For use in High Performance Liquid Chromatography.
Filtered through a 0.2 micron filter.Packaged under Nitrogen.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000