
3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Poly(ethylene glycol) – average mol wt 200

Product Number: **P3015**
CAS Number: 25322-68-3
MDL: MFCD00081839
Formula: (C2H4O)nH2O



TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless to Faint Yellow
Appearance (Form)	Liquid
Brookfield Viscosity	43.0 - 57.0 cps
pH	4.5 - 7.5
5%, At 25 Degrees Celsius	
Water (by Karl Fischer)	≤ 0.5 %
UV Absorbance 215nm	≤ 1.0
50% w/v	
Average Molecular Weight	190 - 210
Recommended Retest Period	----------
2 Years	

Specification: PRD.4.ZQ5.10000009614

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.