Acetonitrile

Baker Analyzed LC-MS Reagent

Material No: 9821.1000GL
Effective from: 2019/07/18

Product Specification

Test	Specification
Assay (by GC)	>= 99.9 %
Residue after Evaporation	<= 1 ppm
Water (H2O)	<= 0.01 %
LC-Gradient-Diode Array Detection (AU), test solution is modified with 0.1% (v/v) formic acid:	
at 220 nm	<= 0.002
at 254 nm	<= 0.001
LC-MS Gradient Suitability Test (TIC, 100 to 2000 m/z), test solution is modified with 0.1% (v/v) formic acid:	
Positive ESI-MS Sensitive Impurities (as Reserpine)	<= 50 ng/ml
Trace Impurities (in ppb)	
Aluminium (Al)	<= 50
Calcium (Ca)	<= 50
Iron (Fe)	<= 50
Magnesium (Mg)	<= 50
Potassium (K)	<= 50
Sodium (Na)	<= 50
Ultraviolet Absorbance (1.00-cm path vs water):	
at 200 nm	<= 0.05
at 220 nm	<= 0.01
at 254 nm	<= 0.01

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000

Material No: 9821.1000GL

For Laboratory, Research or Manufacturing use.
For Use in Liquid Chromatography and Mass Spectrometry
Element concentrations are at time of lot release.



Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000