Specification

Silica Gel 2,5-6 mm with indicator (without cobalt chloride) for analysis, ACS 135571

Physical Description:	Solid
Product Code:	135571
Product Name:	Silica Gel 2,5-6 mm with indicator (without cobalt chloride) for analysis, ACS
Short Description:	additional product description: with moisture indicator
Specifications:	Adsorption capacity (24 h, r.H. = 80 %): min. 27 % Loss on drying (150°C): max. 2.0 % Amber colour: activated Light yellow colour: worn out (regenerate)
WGK:	nwg
Storage:	RT
CAS:	112926-00-8
CS:	28112200

AppliChem GmbH
Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • **info.de@itwreagents.com** • **www.itwreagents.com**
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340