BAKERBOND spe, PCB-A

BAKERBOND spe

Disposable Extraction Columns

Material No: 7511-04
Effective from: 2024/05/28

Product Specification

Test	Specification
Aromatic Sulfonic Acid (SO3H), Strong Cation Exchanger:	
Average Particle Diameter, (APD)	47 - 60 um
Carbon (C)	For Information Only %
Sulfur (S)	For Information Only umol/l
Surface Coverage, Based on %C	180 - 245 ug/m2
Surface Coverage (ug/m2), Based on %S	50 - 80 ug/m2
Extraction Residue	<= 0.2 %
Silica Gel/ Sulfuric Acid (SiOH/H2SO4):	
Identification of Silica treated with Sulfuric Acid	passes test
pH (1 g/100 ml)	<= 2

For Laboratory, Research or Manufacturing use.
3mL Solid Phase Extraction Columns
1000 mg per Column

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000