COVERIS
HIGH PERFORMANCE PACKING

PRODUCT SPECIFICATION
ABSORBENT PADS FOR MEDICAL APPLICATIONS

Customer:	VWR
Product:	Absorbent Pad 50mm x 50mm, 3000cc/m2
Product Code:	FRY-301
Customer Code:	COVSFRY-301
Structure:	Mono-layer structure as follows:
Absorbent Layer:	Proprietary blend of absorbent and bonding fibres, including superabsorbent fibres.
Absorbency:	Nominal absorbency – 3000 cc/m^2
	Nominal absorbency per pad – 7.5 cc/m^3
	Minimum absorbency tolerance is 25% below the nominal
	Absorbency based on test procedures (details on request)
	*Due to processing of material and lamination properties the absorbency of products will on occasion be below the nominal figure. Testing should be conducted to ensure that the pad and required absorbency is suitable for the required application.
Dimensions:	80mm x 50mm
Thickness:	N/A – only absorbency is specified
Weight:	Nominal individual pad weight: 0.2g
Packaging:	Boxed in blue, food-grade liners. 7500 pads per box / 60 boxes per pallet
Palletising:	Boxes shrink-wrapped on standard 1200mm x 1000mm pallet
Expiry Date:	12 Months and must be stored within an ambient temperature.

Absorbent pads are manufactured in ISO 13485:2016 & LSAS accredited facilities
Certification is available upon request.

Disclaimer: This specification is intended as a general guide only. It is recommended that full trials are carried out to ensure the suitability of absorbent pads for individual applications. For confidentiality reasons some of the specification details remain undisclosed.

Coveris Hartlepool
(trading name of Coveris Flexibles (Gainsborough) UK Ltd)
Skerne Road, Oakesway Business Park, Hartlepool TS24 0RH
Tel: +44 (0)1429 272 102 Fax: +44 (0)1429 860 388 Web: www.coveris.com

A member of Coveris Group

Registered in England and Wales No.04786458