epredia Enhancing precision cancer diagnostics

Technical Data Sheet

Product Identification	**Adhesion Microscope slides for in-vitro-diagnostic use**

Article code 1-9							
No. Area	**1**	**2-5**	**3-4**	**5-12**	**13**	**14**	**15**
Art. code digits	**J**	**XXXX**	**A**	**XXXXXXXX**	**X**	**X**	**X**
Description	Adhesion Microscope Slide types	Microscope Slides	Dimension Tolerance	Dimension/Color* *Color for Superfrost microscope slides: white, blue, pink, yellow, green, orange, purple	Frosting	Thickness	Glass Type

Usage	ready to use, single use only
Raw material	extra white glass
Refractive Index	n (λ = 546,07 mm): 1,5171
Abbe value	n/a

Dimensions	length: 75,36 mm (+/- 0,43 mm) width: 25,10 mm (+0 / -0,46 mm)
Thickness	**1,0 mm** (+/- 0,05 mm)

Intended Use & Quality Assurance	
Standards and directive:	See Declaration of Conformity
Quality standards:	1. cleanliness 2. dimensions 3. damages 4. adhesion / stuck 5. packaging
Storage and usage instructions:	See Instruction for use

This document is valid without an original signature. Please contact Epredia Gerhard Menzel GmbH (Quality Management) in any case of questions.

Epredia
Gerhard Menzel GmbH
Saarbrückener Strasse 248 | 38116 Braunschweig
Tel: +49 531 59008-0