FEATURES:

Easy-to-read high contrast graduations

Color-coded magnifier stripes make volume reading and size selection easier (except China manufactured pipets)

Negative graduations

Reverse graduations (except 1mL and 2mL)

Polyester fiber plug

Individually wrapped and bulk packaging options

Extractables – Data is available upon request (with an active NDA in place)

ASPIRATING PIPETS:

Individually wrapped packaging

Extractables – Data is available upon request (with an active NDA in place).

CLAIMS:

Refer to product lot certificate for specific claims information.

LOT NUMBER DESIGNATION:

8 Digit Lot Number: First 3 digits – Julian date, start of manufacturing; Next 2 digits – Year of manufacture; Last 3 digits – Batch identification.

OR

7 Digit Lot Number: First digit – Last number of year of manufacture, Next 3 digits – Julian date, start of manufacturing, Last 3 digits – Batch identification

SHELF LIFE:

5-years from date of manufacture

REGULATORY COMPLIANCE:

Quality Management System – Complies with the current version of the EN ISO 13485 Standard and the FDA CFR 21 PART 820, current Good Manufacturing Practices (cGMP).

This document contains proprietary and confidential information and such information must not be disclosed to others for any purpose or used for manufacturing purposes without written permission from Corning Incorporated.

MANUFACTURING LOCATION DETERMINATION:

Product manufacturing location is available on the Certificate of Compliance.

RESPONSIBLE SITE OF MANUFACTURING	COUNTRY OF ORIGIN
Kennebunk, Maine	United States
Durham, North Carolina	United States
Wujiang, China	China

TABLE I: MATERIALS

COMPONENT (made in USA)	MATERIAL
PIPET	POLYSTYRENE (PS)
PLUG	POLYESTER FIBER
INK	BLACK GRADUATION, COLOR BAND

COMPONENT (made in China)	MATERIAL
PIPET	POLYSTYRENE (PS)
PLUG	POLYETHELYNE TEREPHTHALATE (PET)
INK	BLACK GRADUATION, COLOR BAND

TABLE II: COLOR INK ASSIGNMENT

SIZE (mL)	US MFG
1	YELLOW
2	GREEN
5	BLUE
10	ORANGE
25	RED
50	PURPLE
100	AQUA

TABLE III: PRODUCT AND PACKAGING CONFIGURATION

SIZE (mL)	INCREMENTS	OVERFILL CAPACITY	REVERSE GRAD	CAT. NO.	PACKAGING	CASE QTY	PACK QTY
1	1/100 mL	0.2 mL		4010	BULK	1000	25
				4011	PLASTIC	1000	50
				4012	PLASTIC	200	50
				4485	PAPER-PLASTIC	1000	50
				7041	PAPER-PLASTIC TRIPLE BAG	1000	50
2	1/100 mL	0.8 mL		4020	BULK	1000	25

This document contains proprietary and confidential information and such information must not be disclosed to others for any purpose or used for manufacturing purposes without written permission from Corning Incorporated.

				4021	PLASTIC	1000	50
				4486	PAPER-PLASTIC	1000	50
				7042	PAPER-PLASTIC TRIPLE BAG	1000	50
5	1/10 mL	2.5 mL	√	4050	BULK	500	50
				4051	PLASTIC	200	50
				4487	PAPER-PLASTIC	200	50
				7011	PLASTIC TRIPLE BAG (W/O PLUG)	200	50
				7045	PAPER-PLASTIC TRIPLE BAG	200	50
10	1/10 mL	3.0 mL	√	4100	BULK	50	500
				4101	PLASTIC	50	200
				4492	PLASTIC WIDE TIP	50	200
				4488	PAPER-PLASTIC	50	200
				7015	PAPER-PLASTIC TRIPLE BAG	50	200
				7012	PLASTIC TRIPLE BAG (W/O PLUG)	50	200
25	2/10mL	10.0 mL	√	4250	BULK	25	200
				4251	PLASTIC	50	200
				4489	PAPER-PLASTIC	25	200
				7016	PAPER-PLASTIC TRIPLE BAG	25	200
				7013	PLASTIC TRIPLE BAG (W/O PLUG)	25	200
				7044	PLASTIC TRIPLE BAG (W/O PLUG) LOW PARTICULATE TESTED	25	200
50	1/2mL	10.0 mL	√	4500	BULK	25	100
				4501	PLASTIC	25	100
				4490	PAPER-PLASTIC	25	100
				7017	PAPER-PLASTIC TRIPLE BAG	25	100
				7015	PLASTIC TRIPLE BAG (W/O PLUG)	50	200
				7046	PLASTIC TRIPLE BAG (W/O PLUG) LOW PARTICULATE TESTED	10	100

This document contains proprietary and confidential information and such information must not be disclosed to others for any purpose or used for manufacturing purposes without written permission from Corning Incorporated.

100	1mL	N/A	√	4484	PLASTIC	10	100
				4491	PAPER-PLASTIC	10	100
				7000	PAPER-PLASTIC TRIPLE BAG	10	100
				7500	PLASTIC TRIPLE BAG LOW PARTICULATE TESTED	10	100
				7443	WIDE TIP	10	100

ASPIRATING PIPETS:

SIZE (mL)	INCREMENTS	OVERFILL CAPACITY	REVERSE GRAD	COLOR	CAT. NO.	PACKAGING	CASE QTY	PACK QTY
1	N/A	N/A	N/A	N/A	4975	PAPER-PLASTIC	1000	25
2	N/A	N/A	N/A	N/A	9016	PAPER-PLASTIC	500	50
2	N/A	N/A	N/A	N/A	9186	PLASTIC	1000	100
5	N/A	N/A	N/A	N/A	9099	PLASTIC	200	50

TABLE IV: PRODUCT DIMENSIONS

US MANUFACTURED

Size (mL)	A	B	C	D	E	F	Plant
1	279.40 11.00"	NA	4.65 0.183"	4.65 0.183"	1.68 0.066"	2.97 0.117"	Durham
2	280.67 11.05"	NA	6.27 0.247"	6.27 0.247"	1.63 0.064"	4.17 0.164"	Durham
5	346.71 13.65"	2.57 0.101"	8.00 0.315"	8.00 0.315"	1.93 0.076"	5.97 0.235"	Kennebunk
10	344.30 13.56"	3.51 0.138"	7.75 0.305"	9.53 0.375"	1.60 0.063"	5.46 0.215"	Kennebunk
10 Wide Tip	332.71 13.10"	5.31 0.209"	7.75 0.305"	9.53 0.375"	3.51 0.138"	5.46 0.215"	Kennebunk
SKU 4488	344.93 13.58"	3.30 0.130"	7.62 0.300"	9.60 0.378"	1.52 0.060"	5.72 0.225"	Durham
25	343.92 13.54"	4.70 0.185"	7.75 0.305"	14.35 0.565"	3.18 0.125"	5.46 0.215"	Kennebunk
50	368.30 14.50"	5.49 0.216"	7.62 0.300"	17.78 0.700"	3.18 0.125"	5.08 0.200"	Durham
100	345.74 13.61"	10.46 0.412"	7.75 0.305"	24.38 0.960"	7.14 0.281"	5.46 0.215"	Kennebunk

This document contains proprietary and confidential information and such information must not be disclosed to others for any purpose or used for manufacturing purposes without written permission from Corning Incorporated.

Legend

A	B	C	D	E	F
Length	Tip OD	Mouthpiece OD	Pipet OD	Tip ID	Mouthpiece ID

CHINA MANUFACTURED


Size (mL)	A	B	C	D	E	F
1	274	2.4	4.7	4.7	1.35	2.95
	10.787”	0.094”	0.185”	0.185”	0.053”	0.116”
2	280	2.5	5.74	5.74	1.5	3.92
	11.024”	0.098”	0.226”	0.226”	0.059”	0.154”
5	344	2.7	8	8	1.7	6
	13.543”	0.106”	0.315”	0.315”	0.067”	0.236”
10	341	3.4	7.8	10	1.55	5.25
	13.425”	0.134”	0.307”	0.394”	0.061”	0.207”
25	308.81	4.78	7.62	15.16	2.36	5.08
	12.158”	0.188”	0.300”	0.597”	0.093”	0.200”
50	369	5.35	8.2	17.8	2.9	5.65
	14.528”	0.211”	0.323”	0.701”	0.114”	0.222”
100	345.74 13.61"	10.46 0.412"	7.75 0.305"	24.38 0.960"	7.14 0.281"	5.46 0.215"

Legend

A	B	C	D	E	F
Length	Tip OD	Mouthpiece OD	Pipet OD	Tip ID	Mouthpiece ID

This document contains proprietary and confidential information and such information must not be disclosed to others for any purpose or used for manufacturing purposes without written permission from Corning Incorporated.

REPRESENTATIVE DIAGRAMS:



Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic process.

This document contains proprietary and confidential information and such information must not be disclosed to others for any purpose or used for manufacturing purposes without written permission from Corning Incorporated.

Revision History				
Rev.	Revision Date	Description of Change	Change No.	Revised By
1	08/07/2026	Original Issue	CC-30933	Corinne Yartym

This document contains proprietary and confidential information and such information must not be disclosed to others for any purpose or used for manufacturing purposes without written permission from Corning Incorporated.