Identification of the product

Product number	4901003
Product type	Glass beads
GMDN - Code	none
GMDN - Definition	none
HS-Code	70021000
UDI-DI (smallest sales unit)	04250317312334
UDI-DI (Multipack)	04250317393968

Manufacturer / Supplier Identification

Trademark owner	Paul Marienfeld GmbH & Co.KG Am Wöllerspfad 4, 97922 Lauda-Königshofen, Germany
Country of origin	Czech Republic

Specification & Compliance

Norm	
Diameter	3 mm Ø
Material	soda lime glass
Characteristics	solid
Packing unit per box	500 g in a polybag
Mutipack	20 polybags of 500 g in a carton
Brand	MARIENFELD SUPERIOR

Packaging

Smallest sales unit	0.5 pieces
Multipack	10 pieces