Cole-Parmer®

Cole-Parmer® Gelation Timers - GT-200-5, GT-200-6

- Digital timer: 2 models available, 1rpm and 10rpm
- Stops automatically when gelation occurs
- Easy to use, push button display
- Can be mounted on a retort stand for ease of use
- Plungers are removable for easy cleaning
- Sample cups and disposable plungers are available
- 3 year warranty



Gelation Timer GT-200-5

coleparmer.com

antylia SCIENTIFIC



Gelation Timer GT-200-6

Cole-Parmer® Gelation Timers - GT-200-5, GT-200-6

With the Cole-Parmer® Gelation Timer, gelation time is the transition from liquid to solid during polymerization. Test requires approximately 100 mL of sample Gelation timers are used to measure the transition from liquid to solid during polymerisation. The Techne® Gelation Timer is the ideal instrument for accurate measurement and quality control for the gelation of resin or adhesive based samples in the laboratory.

All models have digital timers, with the time in 1 minute increments (GT-5) up to 9999 minutes or in tenths of a minute (GT-6) up to 999.9 minutes. Every unit is supplied with a calibration certificate for traceability. A bright LED display shows the time to the nearest minute or tenth of a minute, and an audible bleep and light confirm completion of gelation.

Technical Specification

Specification	GT-200-5	GT-200-6
Digital timer increment	1 minute	1/10th minute
Maximum time	9999 minutes	999.9 minutes
Dimensions (L x W x H), mm	94 x 69 x 119	94 x 69 x 119
Voltage	230V, 50Hz	230V, 50Hz
Power, W	5W	5W
Shipping Weight, kg	1.32	1.32

Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Cole-Parmer® GT-200-5 gelation timer, 1opm (230V, 50Hz) & steel plunger	98930-43	GT-200	GT-200-5	FGT5
Cole-Parmer® GT-200-5-120 gelation timer, 1opm (120V, 60Hz) & steel plunger	98930-42	GT-200	GT-200-5-120	FGT5 120V 60
Cole-Parmer® GT-200-6 gelation timer, 10opm (230V, 50Hz) & steel plunger	98930-41	GT-200	GT-200-6	FGT6
Cole-Parmer® GT-200-6-120 gelation timer, 10opm (120V, 60Hz) & steel plunger	98930-40	GT-200	GT-200-6-120	FGT6 120V 60

USA:
+1.800.323.4340
+1.847.549.7600

Canada: +1.800.363.5900
China: +86.21.5109.9909
France: +33 (0) 1486 37800
Germany: +49 (0) 9377 92030

India: +1.800.266.1244
Italy: +39 (0) 1313 89513
UK: +44 (0) 1480 272279
All other countries: +1.847.549.7600

Cole-Parmer® essentials
an antylia scientific company

AZ-EMEA-LT-00299 0922

Find out more at antylia.com