Technical Datasheet

# VWR WIPES, CLEANSMOOTH, WHITE

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

*VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document*

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

**Applicable part numbers :** 115-2394, 115-2395

| Attribute | Value* | Test method |
|---|---|---|
| Composition : | | |
| wood pulp | 55% | |
| polyester | 45% | |
| Production method | spunlace | |
| Colour | white | |
| Basis weight | 68 g/m² | |
| Absorption | 400% | |
| Tensile strength : | | |
| Machine direction | 155,7 N | |
| Cross direction | 85,0 N | |

**Characteristics :**
Low fiber and particle release ; Very strong, even wet ; Resistant to acids, solvents and chlorine ; Good absorption. Recyclable. Food Contact Certified. Does not contain binders. Silicone-free.

**Applications :**
Pharma ; Aerospace ; Semiconductor ; Food & Beverage, Automotive, Laboratories ; Printing, Healthcare.

**Formats & Packaging :**

| Cat No | Type | Format size | Quantity | Primary pkg | Secondary pkg | Pallet |
|---|---|---|---|---|---|---|
| 115-2394 | Industrial roll | 38x30cm | 2 rolls x 500 pcs | Plastic film | n/a | 42 x 2rolls |
| 115-2395 | Z-folded wipes | 38x30cm | 10 packs x 50 pcs | Plastic bag | Cardboard case | 64 cases |

**** Disclaimer:*** *The information presented here is the result of testing and other evaluations conducted during the development of the product. This information does not constitute a specification.*

Date of issue: **March 2017**